


Product Datasheet

G protein-coupled receptor family C group 6 member A precursor antibody (orb73565)

Description	G protein-coupled receptor family C group 6 member A precursor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q5T6X5
Expiration Date	12 months from date of receipt.



Western
Blot
analysis
of A:Raji
whole
ce...