



Product Datasheet

Cofilin 2 antibody (orb73585)


Description	Cofilin 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y281
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Raji whole cell lysate.



Immunohistochemical analysis of formalin-fixed tissue.



Confocal immunofluorescence analysis of cells.