



Product Datasheet

beta IV Tubulin antibody (orb73624)

Description	beta IV Tubulin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P04350
Expiration Date	12 months from date of receipt.



Western Blot
analysis of A:293T
whole ce...



Confocal
immunofluorescence
analysis of ...