



Product Datasheet

DNAJB6 antibody (orb73673)

Description	DNAJB6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	075190
Expiration Date	12 months from date of receipt.



Western Blot
analysis of A:Hela
S3, B:Mo...



Confocal
immunofluorescence
analysis of ...