

## Product Datasheet

### EPCAM antibody (orb736759)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to EPCAM |
| <b>Species/Host</b>        | Mouse                              |
| <b>Reactivity</b>          | Human                              |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | IHC                                |
| <b>Target</b>              | EPCAM                              |
| <b>Note</b>                | For research use only              |
| <b>Isotype</b>             | IgG1                               |
| <b>Clonality</b>           | Monoclonal                         |
| <b>MW</b>                  | 34932                              |
| <b>Uniprot ID</b>          | <b>P16422</b>                      |
| <b>Entrez</b>              | <b>4072</b>                        |
| <b>Dilution Range</b>      | IHC: Ready-to-use                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |