

Product Datasheet

SARS-CoV-2 S1 RBD ELISA Kit (orb736882)

Description	Novel Coronavirus Spike Glycoprotein RBD (SARS-CoV-2 S1 RBD)
Range	1.25 ng/mL-80 ng/mL
Target	Novel Coronavirus Spike Glycoprotein RBD (SARS-CoV-2 S1 RBD)
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, swabs
Uniprot ID	P0DTC2
Sensitivity	0.31 ng/mL
Expiration Date	Please enquire.