

Product Datasheet

MKI67 antibodpy (orb736893)

Description	Rabbit polyclonal antibody to MKI67
Species/Host	Rabbit
Reactivity	Hamster
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	MKI67
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	358694
Uniprot ID	P46013
Dilution Range	WB: 1:500, WB: 1:1000
Expiration Date	12 months from date of receipt.