

Product Datasheet

H1-5 antibodpy (orb736918)

Description	Rabbit polyclonal antibody to H1-5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P16401
Dilution Range	WB: 1/500 - 1/1000, IH: 1/50 - 1/200, IF: 1/50 - 1/200
Expiration Date	12 months from date of receipt.