

Product Datasheet

MARK1 antibopdy (orb737145)

Description	Rabbit polyclonal antibody to MARK1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9P0L2
Dilution Range	WB: 1/500 - 1/1000, IH: 1/50 - 1/200, IF: 1/50 - 1/200
Expiration Date	12 months from date of receipt.