



Product Datasheet

20S Proteasome alpha3 antibody (orb73724)

Description	20S Proteasome alpha3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P25788
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Molt-4 whole ...



Immunohistochemical analysis of formalin...