

## Product Datasheet

### DEGP1 antibopdy (orb737400)

|                                |                                        |
|--------------------------------|----------------------------------------|
| <b>Description</b>             | Rabbit polyclonal antibody to<br>DEGP1 |
| <b>Species/Host</b>            | Rabbit                                 |
| <b>Reactivity</b>              | A. thaliana                            |
| <b>Conjugation</b>             | Unconjugated                           |
| <b>Tested<br/>Applications</b> | WB                                     |
| <b>Note</b>                    | For research use only                  |
| <b>Clonality</b>               | Polyclonal                             |
| <b>Uniprot ID</b>              | <b>O22609</b>                          |
| <b>Dilution Range</b>          | WB: 1/500 - 1/2000                     |
| <b>Expiration Date</b>         | 12 months from date of receipt.        |