

Product Datasheet

Human DEFB4 protein (orb737421)

Description	Recombinant human beta-Defensin 4 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 um filtered solution of 20mM PB, 130mM NaCl, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 130mM NaCl, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Glu23-Pro72
Source	E.coli
Uniprot ID	Q8WTQ1
Expiration Date	6 months from date of receipt.