

Product Datasheet

Human UPK2 protein (orb737427)

Description	Recombinant human Uroplakin-2 protein
Endotoxins	Less than 0.1 EU/ug
Preservatives	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Protein Sequence	Asp26-Gly155
Source	Human Cells
Uniprot ID	O00526
Expiration Date	6 months from date of receipt.