

Product Datasheet

Human IL-3R alpha protein (orb737431)

Description	Recombinant human Interleukin-3 Receptor Subunit Alpha protein
Endotoxins	Less than 0.1 EU/ug
Preservatives	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-hFc
Note	For research use only
Protein Sequence	Lys20-Arg305
Source	Human 293 Cells
Uniprot ID	P26951
Expiration Date	6 months from date of receipt.