



Product Datasheet

SQSTM1 antibody (orb73744)

Description	SQSTM1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13501
Expiration Date	12 months from date of receipt.



Western Blot
analysis of A:293T,
B:Raji ...



Confocal
immunofluorescence
analysis of ...