


## Product Datasheet

### voltage-dependent calcium channel gamma-5 subunit isoform b antibody (orb73746)

|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | voltage-dependent calcium channel gamma-5 subunit isoform b antibody                      |
| <b>Species/Host</b>        | Rabbit                                                                                    |
| <b>Reactivity</b>          | Human                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | WB                                                                                        |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                     |
| <b>Clonality</b>           | Recombinant                                                                               |
| <b>Uniprot ID</b>          | <b>Q9UF02</b>                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |



Western  
Blot  
analysis of  
A:HepG2  
whole c...