

Product Datasheet

Human CELA3A protein (orb737469)

Description	Recombinant human Chymotrypsin-Like Elastase Family Member 3A protein
Endotoxins	Less than 0.1 EU/ug
Preservatives	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10%Glycerol, pH 7.5
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10%Glycerol, pH 7.5
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Protein Sequence	Ser16-His270
Source	Human 293 Cells
Uniprot ID	P09093
Expiration Date	6 months from date of receipt.