

## Product Datasheet

### POU2F2 Immunizing peptide (orb737717)

|                         |                                                                     |
|-------------------------|---------------------------------------------------------------------|
| <b>Description</b>      | POU2F2 Immunizing peptide                                           |
| <b>Preservatives</b>    | 100µg of Lyophilised peptide,<br>formulation De-ionized (DI) water. |
| <b>Storage</b>          | Shipped at ambient temperature, store<br>at -20°C                   |
| <b>Note</b>             | For research use only                                               |
| <b>Protein Sequence</b> | SPTGPSTKIKAEED                                                      |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                     |