

## Product Datasheet

### CENPE Immunizing peptide (orb738106)

|                         |                                                                     |
|-------------------------|---------------------------------------------------------------------|
| <b>Description</b>      | CENPE Immunizing peptide                                            |
| <b>Preservatives</b>    | 100µg of Lyophilised peptide,<br>formulation De-ionized (DI) water. |
| <b>Storage</b>          | Shipped at ambient temperature, store<br>at -20°C                   |
| <b>Note</b>             | For research use only                                               |
| <b>Protein Sequence</b> | KREAHKQVTCEN                                                        |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                     |