


Product Datasheet

CD41/Integrin alpha 2b antibody (orb73821)

Description	CD41/Integrin alpha 2b antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08514
Expiration Date	12 months from date of receipt.



Western
Blot analysis
of A:A549
whole ce...