



Product Datasheet

PSR antibody (orb73822)

Description	PSR antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q6NYC1
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431 whole ce...



Immunohistochemical analysis of paraffin...