

Product Datasheet

BAK/BAK1 Antibody (orb738368)

Description

BAK/BAK1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ICC, IF, IHC, WB

Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human BAK1 (74-95aa QVGRQLALIGDDINRRYDSEFQ), different from the related mouse and rat sequences by one amino acid.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

23409 MW

Uniprot ID


[Q16611](#)

Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, Rat
Immunocytochemistry/Immunofluorescence , 5µg/ml, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of BAK/BAK1 using ...



IHC analysis of BAK/BAK1 using anti-BAK/...



IHC analysis of BAK/BAK1 using anti-BAK/...