

Product Datasheet

NEFL/NF-L Antibody (orb738390)

Description

NEFL/NF-L Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications

Immunogen

E.coli-derived human NEFL/NF-L recombinant protein (Position: F4-A463).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

72 kDa

Uniprot ID


[P07196](#)

Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Mouse, Rat Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Mouse, Rat Immunofluorescence, 5 µg/ml, Mouse, Rat Flow Cytometry, 1-3µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC analysis of NEFL/NF-L using anti-NEF...



Flow Cytometry analysis of 293T cells using...



Western blot analysis of NEFL/NF-L using...