

Product Datasheet

Transketolase/TKT Antibody (monoclonal, 2I3) (orb738409)

Description

Transketolase/TKT Antibody (monoclonal, 2I3)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human Transketolase/TKT recombinant protein (Position: M1-A116).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

68 kDa

Uniprot ID


[P29401](#)

Dilution Range


Western blot, 0.1-0.25µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 1-2µg/ml, Human
Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of Transketolase/TKT



IHC analysis of Transketolase/TKT using ...



IHC analysis of Transketolase/TKT using ...