

Product Datasheet

MAG Antibody (monoclonal, 2G11) (orb738410)

Description

MAG Antibody (monoclonal, 2G11)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC, WB

Immunogen

E.coli-derived human MAG recombinant protein (Position: E34-R605).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW

100 kDa

Uniprot ID


P20916

Dilution Range

Western blot, 0.1-0.25µg/ml, Mouse, Rat Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Mouse,Rat Immunofluorescence, 5 µg/ml, Rat Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of MAG using anti-MAG antibody.



IHC analysis of MAG using anti-MAG antibody.



IHC analysis of MAG using anti-MAG antibody.