

Product Datasheet

Tubulin alpha Antibody (monoclonal, 7B12) (orb738419)

Description

Tubulin alpha Antibody (monoclonal, 7B12)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human Tubulin alpha recombinant protein (Position: N18-A403).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW


56 kDa

Uniprot ID
[Q71U36](#)
Dilution Range


Western blot, 0.1-0.25µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of Tubulin alpha u...



IHC analysis of Tubulin alpha using anti...



IHC analysis of Tubulin alpha using anti...