

Product Datasheet

SARS-CoV-2 N/ORF9/Nucleocapsid Antibody (orb738426)

Description

SARS-CoV-2 N/ORF9/Nucleocapsid Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

MSDNGPQNQRNAPRITFGGPSDSTGSNQNGERSGARSKQRRPQGLPNNTASW
FTALTQHGKEDLKFPRGQGVPIINTNSSPDDQIGYYRRATRRIRGGDGKMKDLSP
RWYFYLLGTGPEAGLPYGANKDGIWVATEGALNTPKDHIGTRNPANNAIVLQL
PQGTTLPKGFYAEGSRGGSQASSRSSRSRNSSRNSTPGSSRGTSPPARMAGNG
GDAALALLLLDRLNQLESKMSGKQQQQGQTVTKKSAAEASKKPRQKRTATKA
YNVTQAFGRRGPEQTQGNFGDQELIRQGTQDYKHWPQIAQFAPSASAFFGMSRI
GMEVTPSGTWLTYTGAIKLDDKDPNFKDQVILLNKHIDAYKTFPPTPEKKDKKKK
ADETQALPQRQKKQTVTLLPAADLDDFSKQLQQSMSSADSTQA

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

Uniprot ID
P0DTC9
Dilution Range

ELISA, 0.001-0.1µg/ml, Human

Expiration Date

12 months from date of receipt.