

## Product Datasheet

### SARS-CoV-2 NSP8 Antibody (orb738432)

|                            |                                                                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | SARS-CoV-2 NSP8 Antibody                                                                                                                                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                            |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                      |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                             |
| <b>Immunogen</b>           | AIASEFSSLPSYAAFATAQEAYEQAVANGDSEVVLKKLKKSLNVAKSEFDRDAAM<br>QRKLEKMADQAMTQMYKQARSEDKRAKVTSAMQTMLFTMLRKLDNDALNNIIN<br>NARDGCVPLNIIPLTTAAKLMVVIPDYNTYKNTCDGTTFTYASALWEIQVVDAD<br>SKIVQLSEISMDNSPNLAWPLIVTALRANSVAVKLQ                |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                                                                       |
| <b>Concentration</b>       | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.                                                                                                                                                         |
| <b>Storage</b>             | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.                             |
| <b>Note</b>                | For research use only                                                                                                                                                                                                             |
| <b>Application notes</b>   | Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml. |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                        |
| <b>Uniprot ID</b>          | <b>P0DTC1, P0DTD1</b>                                                                                                                                                                                                             |
| <b>Dilution Range</b>      | ELISA, 0.001-0.1µg/ml, Human                                                                                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                   |