

Product Datasheet

Poliovirus Receptor/PVR Antibody (orb738463)

Description

Poliovirus Receptor/PVR Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human Poliovirus Receptor/PVR recombinant protein (Position: D28-E331).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

70-80 kDa

Uniprot ID


P15151

Dilution Range


Western blot, 0.1-0.25µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 1-2µg/ml, Human Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of Poliovirus Receptor/PVR



IHC analysis of Poliovirus Receptor/PVR



IHC analysis of Poliovirus Receptor/PVR