

Product Datasheet

TRIM25/EFP Antibody (orb738470)

Description

TRIM25/EGFP Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, IHC, WB

Applications

Immunogen

E.coli-derived mouse TRIM25/EGFP recombinant protein (Position: M1-H590).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

71 kDa

Uniprot ID


[Q61510](#)

Dilution Range


Western blot, 0.25-0.5µg/ml, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5µg/ml, Mouse
Flow Cytometry, 1-3µg/1x10⁶ cells, Mouse, Rat Direct ELISA, 0.1-0.5µg/ml, Mouse

Expiration Date


12 months from date of receipt.



Western blot analysis of TRIM25/EGFP using anti-TRIM25/EGFP antibody.



IHC analysis of TRIM25/EGFP using anti-TRIM25/EGFP antibody.



IHC analysis of TRIM25/EGFP using anti-TRIM25/EGFP antibody.