

Product Datasheet

Mouse BDNF ELISA Kit (orb738498)

Description	Mouse BDNF ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H129-R247
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA, citrate)
Uniprot ID	P21237
Sensitivity	15 pg/ml
Expiration Date	Please enquire.

