

Product Datasheet

Mouse TNF Alpha/Tumor necrosis factor ELISA Kit (orb738503)

Description Mouse TNF Alpha/Tumor necrosis factor ELISA Kit

Reactivity Mouse

Range 15.6 pg/ml - 1,000 pg/ml (cell culture supernatants), 7.8 pg/ml - 500 pg/ml (mouse serum, plasma)

Tested Applications ELISA

Immunogen Expression system for standard: E.coli; Immunogen sequence: L80-L235

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)


Note For research use only

Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)

Uniprot ID **P06804**

Sensitivity 1 pg/ml

Expiration Date Please enquire.



Mouse TNF
alpha
PicoKine
ELISA Kit
Stand...