

Product Datasheet

Human Lipocalin-2/NGAL ELISA Kit (orb738505)

Description	Human Lipocalin-2/NGAL ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q21-G198
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	culture supernatants, serum, plasma (heparin), saliva and urine
Uniprot ID	P80188
Sensitivity	10 pg/ml
Expiration Date	Please enquire.

