

Product Datasheet

Mouse Lipocalin-2/NGAL ELISA Kit (orb738506)

Description	Mouse Lipocalin-2/NGAL ELISA Kit
Reactivity	Mouse
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q21-N200
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine
Uniprot ID	P11672
Sensitivity	10 pg/ml
Expiration Date	Please enquire.

