



Product Datasheet

ASB9 antibody (orb73853)


Description	ASB9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96DX5
Expiration Date	12 months from date of receipt.



Western Blot analysis
of A:Hela S3, B:He...



Immunohistochemical
analysis of formalin...



Confocal
immunofluorescence
analysis of ...