



Product Datasheet

PUF60 antibody (orb73887)

Description	PUF60 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9UHX1
Expiration Date	12 months from date of receipt.



Western Blot analysis
of A:HepG2,
B:MOLT...



Immunohistochemical
analysis of formalin...



Confocal
immunofluorescence
analysis of ...