



Product Datasheet

Aspartoacylase antibody (orb73903)


Description	Aspartoacylase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P45381
Expiration Date	12 months from date of receipt.



Western Blot analysis
of A:Molt-4, B:Raj...



Immunohistochemical
analysis of formalin...



Confocal
immunofluorescence
analysis of ...