



Product Datasheet

USP15 antibody (orb73979)

Description	USP15 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y4E8
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Molt-4 whole ...



Immunohistochemical analysis of formalin...



Costained with Hoechst 33342

Confocal immunofluorescence analysis of ...