



## Product Datasheet

### Ube2L3 antibody (orb73999)


|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | Ube2L3 antibody                                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                    |
| <b>Reactivity</b>          | Human                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | IF, IHC, WB                                                                               |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                     |
| <b>Clonality</b>           | Recombinant                                                                               |
| <b>Uniprot ID</b>          | <b>P68036</b>                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |



Western Blot analysis of A:Hela, B:HepG2...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...