



Product Datasheet

XLalphas antibody (orb74126)

Description	XLalphas antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P84996
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:HepG2, B:Raji...



Immunohistochemical analysis of formalin...