



Product Datasheet

Histone H2A.Z/H2A.F/Z antibody (orb74132)


Description	Histone H2A.Z/H2A.F/Z antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P0C0S5
Expiration Date	12 months from date of receipt.



Western Blot analysis of A: Jurkat whole cell lysate



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue sections



Confocal immunofluorescence analysis of cells