



Product Datasheet

NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa antibody (orb74165)

Description	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 9, 22kDa antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9Y6M9
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:Hela S3, B:He...



Immunohistochemical analysis of formalin...