

Product Datasheet

succinate-CoA ligase, ADP-forming, beta subunit precursor antibody (orb74212)

Description succinate-CoA ligase, ADP-forming, beta subunit precursor antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB


Storage Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note For research use only


Clonality Recombinant

Uniprot ID **Q9P2R7**

Expiration Date 12 months from date of receipt.



Western Blot analysis of A:HepG2, B:Molt...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...