



Product Datasheet

beta Tubulin 3 antibody (orb74224)

Description	beta Tubulin 3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13509
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western Blot analysis of A:Raji whole ce...