



Product Datasheet

Tubulin beta-8 antibody (orb74302)


Description	Tubulin beta-8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q3ZCM7
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:A431, B:H1299...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...