



Product Datasheet

ICAD Antibody (orb74341)

Description	ICAD Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O54786
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse lung tissue...



Immunohistochemical analysis of formalin-fixed tissue...