



Product Datasheet

Snip3L Antibody (orb74418)

Description	Snip3L Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O60238
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...



Immunohistochemical analysis of formalin...