



Product Datasheet

TACI Antibody (orb74437)

Description	TACI Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O14836
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 (A) and U9...



Immunohistochemical analysis of formalin...