


Product Datasheet

IL-22 Receptor Antibody (orb74455)

Description	IL-22 Receptor Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8N6P7
Expiration Date	12 months from date of receipt.



Western blot analysis of human HepG2 cel...