



Product Datasheet

TIRP Antibody (orb74543)

Description	TIRP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q86XR7
Expiration Date	12 months from date of receipt.



Western blot analysis of human (A), mouse (B), and rat (C) samples.



Immunohistochemical analysis of formalin-fixed paraffin-embedded tissue sections.